

# Electric Cast Cutter With Vacuum Cleaner

BJ1201-III



# Electric Cast Cutter With Vacuum Cleaner

The electric cast cutter with vacuum cleaner is designed to cut and remove casts after fixation and reposition in orthopaedic surgery. It has features of good cutting, high speed and safety.

It consists of a cutting handpiece and a vacuum cleaner.

There is two speed control on the handpiece.

The functional, ergonomic handpiece can fit comfortably in the hand, providing the ideal cutting position. Light weight and low noise.

## Handpiece of cutter

High-speed oscillating frequency: 16,000 cpm

Low-speed oscillating frequency: 12,000 cpm

Swing arc: 1mm

Input voltage: DC 24V

Weight: 900g

Noise: 68dB(A)

Length: 235mm

## Medical vacuum cleaner

Max power: 1400VA

Voltage: 110/220V

Frequency: 50/60Hz

Noise: 58dB(A)

Six-layer filtration system can filter 0.27um fine dust, and the filtration efficiency can reach 99.997%.



Company address: Room 606, No. 338, Gonghexin Road, Shanghai, China  
Factory address: 1-D, 2-A, No. 125, Longpan Road, Jiading District, Shanghai  
Tel: +86-21-6528078 66526503  
Fax: +86-21-66527863  
Post code: 200436  
website: [www.bojin-medical.com](http://www.bojin-medical.com)  
Email: [manager@bojin-medical.com](mailto:manager@bojin-medical.com)